



NOTE:

1. THIS IS A ULTRA-HIGH VACUUM CHAMBER (UHV).
2. WHEN MACHINING VACUUM PARTS, USE OF SILICONE AND SULPHUR-BASED CUTTING FLUIDS IS PROHIBITED. USE ONE OF THE FOLLOWING:
A) CIMCOOL 5 STAR 49
B) TRIM SOL
3. ELECTROPOLISHING IS NEEDED BEFORE WELDING. PRIOR TO ELECTROPOLISHING, THE CHAMBER NEEDS TO GO THROUGH A MULTIPLE STEP CLEANING PROCESS INVOLVING DEGREASING, WASHING AND DRY NITROGEN BLOW DOWN. THE CHAMBER VACUUM SIDE SURFACE ROUGHNESS SHALL BE BETTER THAN 63 MICRON IN RMS AFTER ELECTROPOLISHING.
4. WELD SHALL BE GAS TUNGSTEN ARC (GTAW) OR TUNGSTEN INERT GAS (TIG) ON VACUUM SIDE OF JOINTS.
5. VACUUM CHAMBER SHALL BE LEAK TESTED USING A MASS SPECTROMETER WITH MINIMUM SENSITIVITY FOR HELIUM OF 2×10^{-10} STANDARD CC/SEC PER LEAK METER DIVISION. SUPPLIES:
ALCATEL ASM-110TCL
VEECO NGA 500 5000
VEECO MSA-MS-80 OR MS-18
DUPONT CEC 24-120B
6. CALIBRATION OF THE LEAK DETECTOR SENSITIVITY SHALL BE PERFORMED JUST PRIOR TO TESTING. FINAL TEST WILL CONSIST OF SURROUNDING THE CHAMBER (BAGGING) WITH HELIUM. THE CHAMBER WILL BE REJECTED IF A 2% DEFLECTION IN THE MOST SENSITIVE RANGE OF THE LEAK DETECTOR IS SENSED WITHIN 1 MIN.
7. ALL DIMENSIONS IN [] ARE MILLIMETERS AND ARE FOR REFERENCE ONLY.
8. MACHINE FINISH FOR ALL VACUUM SURFACE BEFORE ELECTROPOLISHING.

REFERENCE SOURCE
1 MDC VACUUM PRODUCTS CORP.
23642 CARROT BOULEVARD
HAYWARD, CA. 94545-1651
(800)-443-8817

PARTS LIST			
ITEM	DESCRIPTION	QTY	UNIT
1	TUBING 2.50 O.D. x .065 WALL	304 SST	1
2	TUBING 4.25 O.D. x .120 WALL	304 SST	2
3	TUBING 4.00 O.D. x .120 WALL	304 SST	1
4	TUBING 5.00 O.D. x .120 WALL	304 SST	3
5	TUBING 3.00 O.D. x .120 WALL	304 SST	1
6	FLANGE 6" NOM. ROTATABLE	304 SST	1
7	FLANGE 4-1/2" NOM. ROTATABLE	304 SST	1
8	FLANGE 8" NOM. ROTATABLE	304 SST	3
9	P4105090909-210001	304 SST	1
10	P9-20 COOLING BLOCK	304 SST	1

SEE B.O.M.